

	<h1 style="margin: 0;">THE GOODELL LABORATORY</h1>	
<p style="text-align: center;">Author</p>	<p style="margin: 0;">K. Karen Lin</p>	<p style="margin: 0;">March 30, 2009</p>
<p style="text-align: center;">Title</p>	<p style="margin: 0;">Intracellular Staining of HSCs (sparse cells) with Carrier Cells</p>	
<p style="text-align: center;">Introduction</p>	<p>This protocol describes an intracellular staining technique, adapted from the protocol developed by Dr. Garrett Heffner (Weissman Lab), and is particularly useful for the following scenarios,</p> <ol style="list-style-type: none"> 1) Small number of purified HSCs (~1000 cells) with carriers 2) Surface marker staining followed by intracellular staining 3) Intracellular staining with phospho-proteins <p>The fluorophores are interchangeable. The standard analysers compatible with this particular protocol are a LSRII with a blue laser, red laser, a violet laser, and an UV-laser (BD Biosciences) as well as an Imagestream with comparable lasers (Amnis).</p>	
<p style="text-align: center;">Materials</p>	<ol style="list-style-type: none"> 1. 16% Paraformaldehyde (PFA) solution (Electron Microscopy Sciences, 15710) 2. 100% Ethanol 3. 1% BSA in PBS 4. 30% Hydrogen peroxide solution (Sigma, H1009) 5. Barcode Dyes with desired concentration <ol style="list-style-type: none"> 1) Alexa 488 (Invitrogen A-20000) 2) Pacific Blue (Invitrogen P-10163) 6. Your antibodies <ol style="list-style-type: none"> 1) Note that to date, ONLY FITC and Pacific Blue conjugated-antbodies (among the fluorophore collection in our lab, including FITC, PE, PE-Cy7, PECy5, APC, and APC-Cy7) are found to survive the “quenching process” that is necessary to eliminate the HRP activity of endogenous hydroperoxidase. To perform any intracellular staining in HSC, pre-sorting cells is highly recommended. 2) Fluorophores for the intracellular staining: Alexa647 7. Tyramide Signal Amplification (TSA) Reagents: TSA Kit is available from Invitrogen with Horseradish Peroxidase (HRP) conjugated goat anti-rabbit IgG antibody (TSA Kit with Alexa Fluor 647 Tyramide, Cat No. T20926). 	
<p style="text-align: center;">Protocol</p>		<p><i>Notes</i></p>
	<ol style="list-style-type: none"> 1. Fix the cells-of-interest: sort cells (or cell cultures) into a fixed volume. Add appropriate 16% PFA into cell suspension to 4% PFA. Tap the tube to mix the cells and incubate at RT for 10 minutes (or on ice for 30 minutes) 	

<p>2.</p>	<p>Permeabilization: Add 100% ICED EtOH directly to the tube to make 90% EtOH solution. Incubate on ice for 15 minutes. Wash with PBS, spin cells down.</p>	<p>If you have been sorting cells for the whole day, you could consider stopping the incubation here, and storing the samples in a -20 freezer. Note that you need to mix the EtOH with the cell suspension THOROUGHLY; otherwise the cells will be frozen and disrupted.</p>
<p>3.</p>	<p>Barcoding: A multiple color barcoding strategy could be employed by using combinations of barcoding dyes (Krutzik and Nolan, 2006). In general, cells are resuspended in 30% EtOH. Barcoding dye should be added in a 1/100 ratio (10ul to an 100ul cell suspension) to the desired concentration. Incubate under room temperature for 10 minutes. Wash with 10X volume of BSA-PBS, spin down the cells.</p>	<p>-BSA-PBS= 1% of BSA in PBS -Barcoding dyes usually come in powder forms. These fluorescent dyes should be resuspended in Methanol, and dried back into powders with a Speed-vac. In our lab, the barcoding dyes are aliquoted at 25ug/tube, and stored in the -20°C freezer. For experiments, the barcoding dye aliquots are resuspended in DMSO for the staining process.</p>
<p>4.</p>	<p>In the mean time, wash and spin down the carrier cells as well. Mix these un-labeled carrier cells (~5x 10⁴ cells) with the target cells.</p>	<p>Carrier cells in our lab are splenocytes that have been fixed and permeablized with 90% EtOH. These carriers can be stored in the EtOH solution at -20°C for at least half an year.</p>
<p>5.</p>	<p>Quenching the residual hydrogen peroxidase with 1% H₂O₂ in 1% BSA-PBS buffer. Incubate at room temperature for 30 minutes. Wash with 10X volume of BSA-PBS buffer, spin down (2200rpm, 6 minutes).</p>	

6.	Staining with the 1 st antibody. Incubate under room temperature for 30 minutes. Wash with 10X volume of BSA-PBS buffer, spin down (2200rpm, 6 minutes).	Primary antibodies should be titered for appropriate dilution.
7.	Staining with HRP-conjugated secondary antibody. Incubate at room temperature for 30 minutes. Wash with 10X volume of BSA-PBS, spin down (2200rpm, 6 minutes)	HRP-conjugated secondary antibodies that have been tested: <ol style="list-style-type: none"> 1) goat anti-rabbit antibody 1:400 (Santa Cruz, sc-2065, Lot G3103) 2) goat anti-rabbit antibody 1:200 (Invitrogen, included in the TSA kit, T20926)
8.	Developing color with Tyramide Signal Amplification (TSA) kit: <ol style="list-style-type: none"> 1) Prepare 100X solution A: add 1ul of 30% H₂O₂ into 200 ul of H₂O₂ 2) Prepare FRESH solution B: -1X solutionA -1X Tyramide (100X dilution) - in amplification buffer (100ul/reaction) 3) Staining with 100ul/reaction of solution B, incubate for 5 minutes, stop reaction by adding 10X volume of BSA-PBS 4) Spin down (2200rpm, 6 minutes), resuspend in 350ul BSA-PBS for flow analysis. 	

References.

Krutzik, P.O., and Nolan, G.P. (2006). Fluorescent cell barcoding in flow cytometry allows high-throughput drug screening and signaling profiling. *Nature methods* 3, 361-368.